

## The Microscope

The microscope led to the discovery of the cell

### Types of Microscopes

#### 1. *Stereomicroscope/Dissecting Microscope*

- a. 2 eyepieces: does not flip the specimen in the field of view
- c. Used to observe fine detail

#### 2. *Compound Light Microscope*

- a. 2 lenses used at a time (eyepiece & objective)
- b. Objects must be very thin for light to pass through

#### 3. *Electron Microscope (most detail)*


- a. Uses electrons to magnify an image.

**Field of View (FOV):** What you see in the illuminated circle

Objects under a microscope are

- Upside down
- Backwards
- Larger/magnified/bigger
- More detailed



**Movement of slide:** In the FOV the specimen moves in the opposite direction of the slide.



Letter	Name of Part	Function
A	Eyepiece/Ocular	Magnifies
B	Coarse Adjustment	Moves stage only when using low power
C	Fine Adjustment	Used to focus

D	Arm	Used to carry the scope
E	Stage Clips	Secure the slide
F	Stage	Slide is placed here
G	Objective	Magnifies
H	Diaphragm	Controls the amount of light
I	Light source	Provides light
J	Base	Used to carry the scope

**Field of view differences under low & high power**

Low Power	High Power
<ul style="list-style-type: none"> <li>• Large</li> <li>• Bright</li> </ul> 	<ul style="list-style-type: none"> <li>• Small</li> <li>• Dark</li> </ul> 

**Total Magnification:** How much the specimen is magnified using the 2 lenses.

*Formula: (Ocular) x (Objective) = TM (total magnification)*

**Microscope Rules**

- Always start with the lowest objective
- Always **center under low power** & then switch to high power
- Never use the coarse adjustment (large knob) when using high power, it may break the lens or the slide.